
SP-220 Gas Leak Checkers

Gas Specification List for Fumigant Version

NO.	Gas name	Display	Conc. Display 1 (ppm)	Display 2 (ppm)	Display 3 (ppm)
1	Phosphine	PH3	0.1	2	4.5
2	Methyl bromide	CH3Br	1	20	100
3	Carbon disulfide	CS2	0.1	-	12
4	Methyl iodide	CH3I	1	10	30
5	Hydrogen cyanide	HCN	1	-	-
6	Sulfuryl fluoride	SO2F2	-	- 1	800
7	Ethylene dibromide	C2H4Br2	1	10	30



Gas sensitivity level & Display

SP-220 automatically changes gas sensitivity level depends on target gas.

Gas sensitivity level can be set on 3 (Type FUM) scales. *Refer to right figures.

◆ How to change the level ◆
Push both 「AIR」+「MODE」
buttons at the same time



Note 1) [-]: It means they have sensitivity, but there is no value.

Note 2) Alarm accuracy, response time etc. are confirmed with calibration gas (PH3)

Note 3) The measuring results of the gas selection are intended for guidance only. For the most accurate measurement is for calibration gas (PH3)



Int. Phone: +61 8 6108 0000 Int. Fax: +61 8 9408 1868 Email: info@gastech.com.au

www.gastech.com.au

HEAD OFFICE 24 Baretta Road Wangara, WA 6065 Phone +61 8 6108 0000 SYDNEY OFFICE 21/25, Narabang Way Belrose NSW 2085 Phone: +61 2 9451 0054

DISTRIBUTED BY: